

01 Eye Specula


Lidsperrer

- Wire 15
- Wire, adjustable 17
- Adjustable, locking 19
- Aspiration, adjustable, locking 25
- Glaucoma 29



Geuder 01 Eye Specula Wire


G-17020 BARRAQUER
Eye Speculum
fenestrated
14 mm blades
overall length 38 mm




G-17025 BARRAQUER
Eye Speculum
for infants
fenestrated
10 mm blades
overall length 40 mm




G-17023 BARRAQUER
Eye Speculum
for premature infants
fenestrated
7 mm blades
overall length 34 mm



G-17026 BARRAQUER
Eye Speculum
15 mm solid blades
overall length 40 mm




G-32136 BARRAQUER
Eye Speculum
fenestrated
14 mm round blades
overall length 38 mm



Geuder® 01 Eye Specula
Wire


GS-32137

BARRAQUER
Eye Speculum
for infants fenestrated
10 mm round blades
overall length 38 mm




GS-32148

BARRAQUER
Eye Speculum
fenestrated
14.5 mm open blades
overall length 40 mm




GS-32149

BARRAQUER/ FRANKFURT MODEL
Eye Speculum
fenestrated, large
14.5 mm blades
overall length 40 mm



GS-17004


H.-R. KOCH
Temporal Approach Eye Speculum
for nasal placement
fenestrated
14 mm open blades
overall length 40 mm



Geuder® 01 Eye Specula
Wire

GS-32129


Temporal Approach Eye Speculum
for nasal placement
fenestrated
15 mm open blades
overall length 40 mm



Geuder® 01 Eye Specula
Wire, adjustable

GS-17010

BARRAQUER - OOSTERHUIS
Eye Speculum
fenestrated, adjustable
14 mm blades
overall length 40 mm



GS-17015


BARRAQUER - OOSTERHUIS
Eye Speculum
solid blades, adjustable
15 mm blades
overall length 40 mm



GS-32138

GS-32139

BARRAQUER - OOSTERHUIS
Eye Speculum
fenestrated round blades
adjustable
14 mm blades, overall length 40 mm
11 mm blades, overall length 30 mm, for infants



Geuder® 01 Eye Specula

Wire, adjustable

⌘-17017


BARRAQUER - OOSTERHUIS - KRATZ

Eye Speculum

adjustable , fenestrated

15 mm open blades

overall length 40 mm



⌘-32119

KRATZ - BARRAQUER

Eye Speculum

adjustable, fenestrated

15 mm round wing blades

overall length 40 mm



⌘-17027

THOMAS NEUHANN

Eye Speculum

adjustable

15 mm solid blades

overall length 40 mm



⌘-17030

SAUER

Eye Speculum

12 mm solid blades

overall length 30 mm



Geuder® 01 Eye Specula

Wire, adjustable

⌘-32134

TAKEDA

Eye Speculum

for infants

9 mm fenestrated blades

with finger grips

and angled interlacing valves

overall length 28 mm



Geuder® 01 Eye Specula

Adjustable, locking

⌘-15945

MELLINGER


Eye Speculum

large

adjustable, self retaining

16.5 mm fenestrated blades

overall length 65 mm



⌘-32130

MELLINGER


Eye Speculum

small, for infants

adjustable, self retaining


12.5 mm fenestrated blades

overall length 55 mm



Geuder® 01 Eye Specula
Adjustable, locking


G-15946
MELLINGER – BEARD
Eye Speculum
large
adjustable, self retaining
15.5 mm solid blades
overall length 65 mm



G-32135
MELLINGER – BEARD
Eye Speculum
small, for infants
adjustable, self retaining
11.5 mm solid blades
overall length 55 mm




G-16050
G-16055
PRETORI
Eye Speculum
13.5 mm fenestrated blades
adjustable, self retaining
overall length 63 mm
left
right




Geuder® 01 Eye Specula
Adjustable, locking


G-15950
WEISS
Eye Speculum
fenestrated
adjustable with locking screw
17.5 mm blades
overall length 85 mm




G-15951
WEISS
Eye Speculum
fenestrated
adjustable with locking screw
17.5 mm open blades
overall length 85 mm



G-15955
WEISS
Eye Speculum
small
adjustable with locking screw
15 mm fenestrated blades
overall length 55 mm



G-15970
WEISS
Eye Speculum
adjustable with locking screw
16.5 mm solid blades
overall length 85 mm



Geuder® 01 Eye Specula
Adjustable, locking

GS-15975

WEISS

Eye Speculum

small

adjustable with locking screw

14 mm solid blades

overall length 55 mm



GS-16010


BANGERTER

Eye Speculum

adjustable, with locking screw

15 mm solid blades

overall length 52 mm



GS-16025

CASTROVIEJO

Eye Speculum

medium

16 mm blades, fenestrated, adjustable

overall length 96 mm



GS-15990

COOK

Eye Speculum

large

16 mm solid blades

adjustable, with locking screw

overall length 50 mm



Geuder® 01 Eye Specula
Adjustable, locking

GS-15995

COOK

Eye Speculum

medium

11 mm solid blades

adjustable, with locking screw

overall length 40 mm



GS-16000

COOK

Eye Speculum

small

8 mm solid blades

adjustable, with locking screw

overall length 35 mm



GS-32128

COOK

Eye Speculum

for infants

5 mm solid blades

adjustable, with locking screw

overall length 30 mm



GS-17000

GS-17001

MURDOCK

Eye Speculum


14 mm fenestrated blades

adjustable, self retaining

overall length 50 mm

left

right



Geuder® **01 Eye Specula**
Adjustable, locking

GS-17002


GS-17003

MURDOCK
Eye Speculum
small
12 mm fenestrated blades
adjustable, self retaining
overall length 45 mm
left
right



GS-15960

LIEBERMAN - TENNANT
Eye Speculum
adjustable, fenestrated
for temporal placement
15 mm open blades
overall length 80 mm




GS-15965

KERSHNER
Eye Speculum
adjustable, reversible
for nasal and temporal placement
15 mm fenestrated blades
overall length 75 mm



GS-15966

KERSHNER
Eye Speculum
adjustable, reversible
for nasal and temporal placement
15 mm solid blades
overall length 75 mm



GS-33965


DICK
Femtocataract Eye Speculum
adjustable
18 mm flattened fenestrated blades
with middle bar to control redundant lid tissue
overall length 80 mm



Geuder® **01 Eye Specula**
Adjustable, locking

GS-33971

SOLOMON
Femtocataract Eye Speculum
14.5 mm flattened fenestrated blades
with middle bar to control redundant lid tissue
double knurled screw for fast adjustment
and fine tuning
overall length 95 mm



Geuder® **01 Eye Specula**
Aspiration, adjustable, locking

GS-33919


Spare Tubing Set
for eye speculum with aspiration




Scale 1:2

GS-33920


KNORZ
LASIK Eye Speculum
adjustable, with aspiration
16 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adapter



Geuder® 01 Eye Specula
Aspiration, adjustable, locking



NODA
Eye Speculum
adjustable, with aspiration
to remove excess fluid from the surgical field
supplied with silicone tubing and adapter
15 mm open blades
12 mm open blades, for infants




PAULMANN
Eye Speculum
adjustable, with aspiration for PPV,
with extra aspiration bar
16 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adapter




Clamp for Infusion Tube
for PAULMANN Eye Speculum ()

Geuder® 01 Eye Specula
Aspiration, adjustable, locking




LIEBERMAN - TENNANT
Eye Speculum
adjustable, with aspiration for temporal placement
15 mm open blades
overall length 80 mm
to remove excess fluid from the surgical field
supplied with silicone tubing and adapter




KOGA - NODA
Eye Speculum
with 5 mm aspiration slits
adjustable
12 mm open blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adapter



Geuder® 01 Eye Specula
Aspiration, adjustable, locking




Ge-17031
NODA
Eye Speculum
adjustable, with aspiration
open blades
2 x 10 aspiration ports
supplied with silicone tubing and adapter




Ge-33925
MAUS
LASIK Eye Speculum
adjustable, with middle aspiration bar
16 mm fenestrated blades
to remove excess fluid from the surgical field
supplied with silicone tubing and adapter

Geuder® 01 Eye Specula
Glaucoma



Ge-15920
Eye Speculum
fenestrated blades
adjustable mechanism regulates opening and blade elevation
optimal for glaucoma surgery
large, 16 mm blades



G-15921 SCHOTT
Eye Speculum
solid blades
adjustable mechanism regulates opening and blade elevation
optimal for glaucoma surgery
large, 15 mm blades