

# Foreign Body Tools

Fremdkörperinstrumente

Instruments 54  
Needles 55  
Spuds 55



Geuder® 05 Foreign Body Tools  
Instruments



G-14550      **Foreign Body Instrument**  
double-ended with sliding handle  
2.2 mm wide curved needle  
1 mm wide PLANGE spud  
overall length 95 mm



G-14540      **Foreign Body Instrument**  
double ended with sliding handle  
2.2 mm wide curved needle  
1.4 mm wide gouge  
overall length 95 mm

Geuder® 05 Foreign Body Tools  
Needles



G-14490      **Foreign Body Needle**  
straight  
4.5 mm long triangular tip  
2.4 mm wide  
overall length 120 mm



G-14500      **Foreign Body Needle**  
curved  
4.5 mm long triangular tip  
2.4 mm wide  
overall length 120 mm

Geuder® 05 Foreign Body Tools  
Spuds



G-14570      **PLANGE Foreign Body Spud**  
curved and rounded  
1 mm wide spud  
overall length 125 mm