

Geuder®

08 Scissors

Scheren

- Corneal 69

Stitch 71

Synechia 71

Trabeculum 72

Capsule 72
- Iris 75

IOL 76

with Needle Holder, Conjunctival 77

Tenotomy 78

Strabismus 79

Utility, Stitch, Ring Handle 80

Enucleation 83

Vitreoretinal 84



Geuder®

08 Scissors

Corneal

G-19710 G-19715	TROUTMAN Corneal Scissors 7 mm blades 9.5 mm from tip to pivot overall length 100 mm right, curved blades, blunt tips left, curved blades, blunt tips	
G-19690 G-19695	CASTROVIEJO Corneal Scissors lower blade 9.5 mm upper blade 9.0 mm overall length 105 mm right, curved blades, blunt tips left, curved blades, blunt tips	
G-19720 G-19725	KATZIN Keratoplasty Scissors 6 mm blades 8 mm from tip to pivot overall length 100 mm right, curved blades, blunt tips left, curved blades, blunt tips	
G-19718	SAUTTER Corneal Scissors light curve, fine blunt tips 10 mm blades 14 mm from tip to pivot overall length 100 mm	
G-19780	CASTROVIEJO Corneal Scissors strong curve, blunt tips 12 mm blades 15 mm from tip to pivot overall length 103 mm	
G-19790	CASTROVIEJO Corneal Scissors angled to side, blunt tips 12 mm blades overall length 110 mm	

Geuder® 08 Scissors
Corneal

G-19600 G-19605	Corneal Scissors light curve 8 mm blades overall length 105 mm sharp tips blunt tips	
G-19607 G-19615	Corneal Scissors strong curve, blunt tips 6 mm blades, overall length 105 mm 11 mm blades, overall length 110 mm	
G-19727	MACKENSEN Keratoplasty Scissors strong curve, blunt tips overall length 85 mm	
G-19585	Corneal Section Scissors angled up, 8 mm blades, blunt tips overall length 100 mm	

Geuder® 08 Scissors
Stitch

G-19570	Micro Stitch Scissors straight, 8 mm blades, sharp tips overall length 105 mm	
G-19580	FINE Stitch Scissors straight, 14 mm blades, sharp tips overall length 110 mm	
G-19590	CASTROVIEJO Anterior Synechia Scissor angled up 9 mm vaulted blades, blunt tips overall length 105 mm	
G-19740	CASTROVIEJO Anterior Synechia Scissor angled up 15 mm vaulted blades, blunt tips overall length 115 mm	

Geuder® 08 Scissors
Trabeculum

	Trabeculum Scissors sharp tips 8.5 mm blades 1.5 mm pivot width overall length 80 mm	
┐-19745	straight	
┐-19747	curved	

Geuder® 08 Scissors
Capsule

	Capsulotomy Scissors angled up, sharp tips 2 mm pivot width overall length 80 mm	
┐-19776	9 mm blades	
┐-19777	11 mm blades	
┐-32100	Capsulotomy Scissors angled, lightly curved tips 13 mm blades 2 mm pivot width overall length 80 mm	
┐-32101	Capsulotomy Scissors angled, strongly curved tips 13 mm blades 2 mm pivot width overall length 80 mm	
┐-32105	Capsulotomy Scissors lightly curved, sharp tips 11 mm blades 1.8 mm pivot width overall length 88 mm	

Geuder® 08 Scissors
Capsule

┐-19620	TUEBINGEN MODEL Capsulotomy Scissors straight, sharp tips 9 mm blades 1.8 mm pivot width overall length 88 mm	
┐-19630	TUEBINGEN MODEL Capsulotomy Scissors curved, sharp tips 9 mm blades 1.8 mm pivot width overall length 88 mm	
┐-19760	VANNAS Capsulotomy Scissors extra delicate, straight 6 mm blades, sharp tips 1.7 mm pivot width overall length 80 mm	
┐-19765	VANNAS Capsulotomy Scissors extra delicate, curved 6 mm blades, sharp tips 1.7 mm pivot width overall length 80 mm	
┐-19770	VANNAS Capsulotomy Scissors extra delicate, angled to side 7 mm blades, sharp tips overall length 80 mm	
┐-19775	VANNAS Capsulotomy Scissors extra delicate, angled up 7 mm blades, sharp tips overall length 80 mm	

Geuder® 08 Scissors
Capsule


GS-32940

H.-R. KOCH
Capsule Scissors
horizontal, 45° angled
2.5 mm blades, 360° rotatable head
22 Gauge/ 0.7 mm straight shaft
overall length 107 mm



GS-32942

BRAUWEILER
Capsule Scissors
vertical, 1.25 mm blades
22 Gauge/ 0.7 mm curved shaft
overall length 100 mm



GS-38246

YAE
Capsule Scissors
interchangeable, titanium, gold-colored
3 mm straight micro blades
23 Gauge/ 0.6 mm curved shaft
for cutting synechia, iridotomies, and to rescue the rhexis during MICS
overall length 42 mm
for use with JAPAN MODEL Handle GS-38240



GS-38240

JAPAN MODEL
Handle
for capsule scissors and forceps
6 mm, titanium
overall length 90 mm



GS-32948

HATTENBACH
Hybrid Scissors
to cut sutures during iris and
dislocated IOL fixation
1.25 mm blades
25 Gauge/ 0.5 mm angled shaft



Geuder® 08 Scissors
Iris

GS-19540
GS-19550

NOYES
Iris Scissors
straight
18 mm blades
overall length 120 mm
sharp tips
blunt tips



GS-19640
GS-19655

WECKER
Iris Scissors
11 mm blades
overall length 110 mm
sharp tips
blunt tips



GS-19660
GS-19675

WECKER
Iris Scissors
7 mm blades
overall length 110 mm
sharp tips
blunt tips



GS-19680

BARRAQUER
Iris Scissors
7 mm blades
blunt tips
overall length 56 mm




Geuder® 08 Scissors
IOL

GS-32945

STEINMETZ
IOL Scissors


6.5 mm upper blade
serrated lower tip with groove
scissors closed: 0.95 mm wide x 0.75 mm high
overall length 92 mm
for cutting acrylic and silicone foldable IOLs



GS-32946

STEINMETZ
IOL Scissors

for small incision technique
6.5 mm upper blade
serrated lower tip with groove
scissors closed: 0.85 mm wide x 0.65 mm high
overall length 92 mm
for cutting acrylic and silicone foldable IOLs






GS-32998

GS-32996

SNYDER – OSHER
IOL Scissors

for cutting acrylic and silicone foldable IOLs
lower notched blade to create IOL traction
18 Gauge/ 1.2 mm shaft
overall length 150 mm
short model, 3.5 mm blades
long model, 4.5 mm blades



Geuder® 08 Scissors
with Needle Holder, Conjunctival

GS-17440

Combined Needle Holder/ Scissors

curved 12 mm blades with
3 mm needle holder and tying platform
used to hold, place, and cut fine suture material overall
length 125 mm



GS-19750

WESTCOTT
Scissors

blunt tips
curved, standard blades
overall length 120 mm



GS-19751

WESTCOTT
Scissors

pointed tips
curved, standard blades
overall length 120 mm



GS-32250

WESTCOTT
Scissors

blunt tips
curved, delicate blades
overall length 100 mm



Geuder® 08 Scissors
Tenotomy

G-19470
G-19475

STEVENS
Tenotomy Scissors
short blades, blunt tips
22 mm from tip to pivot
overall length 105 mm
straight
curved



G-19480
G-19485

STEVENS
Tenotomy Scissors
long blades, blunt tips
28 mm from tip to pivot
overall length 110 mm
straight
curved



Geuder® 08 Scissors
Strabismus

G-19430
G-19435

BONN MODEL
Strabismus Scissors
blunt tips
24 mm from tip to pivot
overall length 90 mm
straight
curved



G-19420
G-19425

Strabismus Scissors
long blade
28 mm from tip to pivot
overall length 105 mm
straight
curved



G-19460
G-19465

Strabismus Scissors
short blade
28 mm from tip to pivot
overall length 110 mm
straight
curved



Utility, Stitch, Ring Handle

G-19350	90 mm	27 mm
G-19355	100 mm	29 mm
G-19360	115 mm	30 mm
G-19365	120 mm	36 mm



G-19370	90 mm	27 mm
G-19375	100 mm	28 mm
G-19380	110 mm	29 mm
G-19385	120 mm	34 mm



G-19390	straight
G-19395	curved



Utility, Stitch, Ring Handle

G-19400	curved
G-19405	straight



G-19410	SPENCER Stitch Scissors 20 mm straight blades with hook to lift and secure sutures before cutting overall length 90 mm
---------	---



G-19505 **Eye Scissors**
angled to side, probe tip
sharp-blunt
overall length 115 mm



Geuder® 08 Scissors
Utility, Stitch, Ring Handle

G-19510

Eye Scissors
angled to side, probe tips
blunt tips
overall length 115 mm



G-31315

Operating Scissors
straight blades
sharp-blunt
overall length 145 mm



Geuder® 08 Scissors
Enucleation

G-19520
G-19525
G-19530

Enucleation Scissors
curved blades, blunt
40 mm from tip to pivot
overall length 135 mm
strong curve
medium curve
light curve



