

Geuder®

09 Forceps

Pinzetten

- Cilia, Foreign Body 87

Conjunctiva 88

Suturing 89

Keratoplasty 91

Suturing, Tying 91

Fragment, Jewelers 94

Fixation, Dressing, Tissue 95

Iris 95

Capsulorhexis 97

DMEK, DSAEK, LASIK 101

Lenticule 102
- Nucleus Prechoppers, Implantation 103

IOL 104

Superior Rectus 107

Scleral Plug, Trocar 109

Vitreoretinal 110

Chalazion, Ptosis 112

Mosquito 113

Towel, Serrefine 114



Geuder®

09 Forceps

Cilia, Foreign Body

G-18720	HENRY Cilia Forceps straight jaws smooth, blunt, squared tips overall length 100 mm	
G-18730	GRADLE Cilia Forceps curved jaws smooth, blunt, squared tips overall length 98 mm	
G-18750	BARRAQUER Cilia Forceps lightly angled 65 mm platform overall length 109 mm	
G-18740	BEER Cilia Forceps 3 mm platform rounded tips 45 mm wide overall length 78 mm	
G-18850	SHAAF Cilia and Foreign Body Forceps grooved tips curved and pointed finger grip handle overall length 96 mm	
G-31620	BARRAQUER Cilia and Suture Forceps grooved tips, 7 mm platform curved and pointed finger grip handle overall length 112 mm	

Geuder® 09 Forceps
Cilia, Foreign Body


GS-32740

HAMBURG MODEL

Cilia and Foreign Body Forceps

tapered fine tips

overall length 90 mm



Geuder® 09 Forceps
Conjunctiva

GS-31230


KREFELD MODEL

Conjunctiva Forceps

standard

1 x 2 teeth

overall length 100 mm



GS-31235


FREIBURG MODEL

Conjunctiva Forceps

delicate

1 x 2 teeth

overall length 90 mm



GS-18810


DRAEGER

Conjunctiva Forceps

12 mm wide

2 atraumatic fixation knobs

overall length 97 mm



Geuder® 09 Forceps
Suturing

CASTROVIEJO

Suturing Forceps

1 x 2 teeth

6 mm tying platform

overall length 100 mm

0.28 mm, teeth 0.12 mm


GS-19001

GS-19005

GS-19010

GS-19015

GS-19020



CASTROVIEJO

Suturing Forceps

new model

1 x 2 teeth

6 mm tying platform

overall length 100 mm

0.28 mm, teeth 0.12 mm


GS-19022

GS-19023

GS-19024

GS-19025

GS-19026



MUNICH MODEL


Suturing Forceps

1 x 2 teeth, 0.35 mm

55 mm tapered tying platform

overall length 100 mm

GS-19008




Geuder® 09 Forceps
Suturing

18965

18966

ZUERICH MODEL
Suturing Forceps
delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm
straight
angled



Geuder® 09 Forceps
Colibri Style

18950

18960

15134


18975


Colibri Forceps
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 75 mm


Micro Colibri Forceps
delicate
1 x 2 teeth, 0.12 mm
5 mm tying platform
overall length 77 mm

Micro Colibri Forceps
extra delicate
1 x 2 teeth, 0.12 mm
4 mm tying platform
overall length 76 mm

BARRAQUER
Suture Removing Forceps
very fine pointed tips
4 mm tying platform
overall length 77 mm









Geuder® 09 Forceps
Keratoplasty

18972

18970

38623


THOMAS NEUHANN
Keratoplasty Forceps
for corneal suturing and fixation
two prongs with 1 x 2 teeth
0.8 mm prong spread
overall length 80 mm

REMKY
Keratoplasty Forceps
two prongs with 1 x 1 teeth, angled
1.7 mm prong spread
for corneal suturing and fixation
overall length 80 mm

REMKY
DMEK Keratoplasty Forceps
modified by SZURMAN
two prongs with 1 x 1 teeth, straight
1.7 mm prong spread
33 mm shafts
to hold graft during preparation
overall length 89 mm








Geuder® 09 Forceps
Suturing, Tying

18940

18945


PAUFIQUE
Universal Suturing Forceps
1 x 2 teeth
6 mm tying platform
overall length 88 mm
standard
delicate





19045

19046


PAUFIQUE
Tying Forceps
with platform
overall length 88 mm
standard, 6 mm platform
delicate, 5 mm platform



Geuder® 09 Forceps
Suturing, Tying

GS-19040	TUEBINGEN MODEL Tying Forceps curved, smooth jaws 6 mm tying platform overall length 99 mm	 
GS-19035	TUEBINGEN MODEL Tying Forceps delicate, smooth jaws 4 mm tying platform overall length 84 mm	 
GS-19034	BONN MODEL Tying Forceps extra delicate, smooth jaws 3 mm tying platform overall length 73 mm	 
GS-19032	Tying Forceps angled, 7.5 mm 6 mm tying platform overall length 99 mm	 
GS-31625	DARDENNE Tying Forceps angled, 6.5 mm 6 mm tying platform round handle with guide pin overall length 102 mm	 

Geuder® 09 Forceps
Suturing, Tying

GS-19038	BANGERTER Tying Forceps 9.5 mm platform longitudinal V-shape groove overall length 100 mm	 
GS-19048	Tying Forceps utility serrated jaws 1 mm wide overall length 106 mm	 
GS-32947	HATTENBACH Hybrid Tying Forceps angled, for the anterior chamber intraocular 1.5 mm micro tying tips 25 Gauge/0.5 mm overall length 98 mm for use in combination with HATTENBACH Hybrid Needle Holder GS-32949	 

Geuder® 09 Forceps
Fragment, Jewelers

G-31623

KELMAN – McPHERSON
Fragment Forceps
angled, smooth jaws
7.5 mm platform
overall length 100 mm

G-33941

KELMAN – McPHERSON
Fragment Forceps
angled, fine grooved jaws,
for Femto-LASIK to retrieve lenticule
during SMILE procedures
7.5 mm platform
overall length 100 mm

G-31622

G-31624

GASKIN
Fragment Forceps
angled
with tying platform
9 mm platform, overall length 107 mm
11 mm platform, overall length 100 mm

G-37700

Jewelers Type Forceps
standard, straight
overall length 120 mm

G-37705

Jewelers Type Forceps
delicate, straight
overall length 115 mm

G-37710

Jewelers Type Forceps
extra delicate, straight
overall length 115 mm

Geuder® 09 Forceps
Fixation, Dressing, Tissue

G-19050

GRAEFE
Fixation Forceps
without lock
4.5 mm wide
8 interlocking teeth
overall length 110 mm

G-18780

G-18781

G-18782

G-18783

G-18784

Dressing Forceps
serrated jaws
overall length 100 mm

G-18790

G-18791

G-18792

G-18793

G-18794

Tissue Forceps
1 x 2 teeth
overall length 100 mm

Geuder® 09 Forceps
Iris

G-18855

G-18856

Iris Forceps
serrated, delicate
overall length 73 mm
straight
curved

G-18908

G-18911

BONN MODEL
Iris Forceps
serrated, extra delicate
overall length 73 mm
straight
angled at 8 mm


Geuder® 09 Forceps
Iris

┐-18900

┐-18905

┐-18910


BONN MODEL
Iris Forceps
1 x 2 teeth, 0.12 mm
extra delicate
overall length 73 mm
straight
angled at 6 mm
angled at 8 mm



┐-18845


┐-18890

BONN MODEL
Iris Forceps
long, straight
1 x 2 teeth, 0.12 mm
overall length 102 mm
extra delicate
delicate, with 45 mm tying platform




┐-18840

SWISS MODEL
Iris Forceps
curved, 1 x 2 teeth
overall length 100 mm




┐-18857

Iris Forceps
straight, delicate
1 x 2 teeth
overall length 72 mm




┐-31953

Iris Forceps
straight, strong
1 x 2 teeth
overall length 72 mm



┐-32240

SWISS MODEL
Iris Forceps
straight, extra delicate
1 x 2 teeth
overall length 70 mm



Geuder® 09 Forceps
Iris


┐-18897

ELZE
Iris Forceps
extra delicate
1 x 2 teeth
overall length 78 mm



┐-18913

THOMAS NEUHANN
Iris Forceps
curved shaft
1 x 2 teeth
23 Gauge/ 0.6 mm
overall length 95 mm



Geuder® 09 Forceps
Capsulorhexis

┐-31307

UTRATA
Capsulorhexis Forceps
strong, with heel stop to
protect the iris
interlocking micro serrated tips
12 mm straight shafts, angled
overall length 103 mm



┐-31308

Capsulorhexis Forceps
round handle with guide pin
extra delicate with
interlocking micro serrated tips
13 mm curved shafts, angled
overall length 105 mm



┐-31300

MASKET
Capsulorhexis Forceps
delicate, with heel stop to
protect the iris
interlocking micro serrated tips
12 mm curved shafts, angled
overall length 103 mm



Geuder® 09 Forceps
Capsulorhexis

GS-31304

Capsulorhexis Forceps
new model, with guide pin
extra delicate with
interlocking micro serrated tips
13 mm curved shafts, angled
overall length 103 mm





GS-31301

Capsulorhexis Forceps
extra delicate
with smooth triangular grasping tips
13 mm curved shafts, angled
overall length 98 mm





GS-31305

Capsulorhexis Forceps
extra delicate with
interlocking micro serrated tips
13 mm curved shafts, angled
overall length 98 mm






GS-31303

H.-R. KOCH
Capsulorhexis Forceps
extra delicate with
interlocking micro serrated tips
13 mm curved shafts
1 mm marks from tip to bent
overall length 98 mm








GS-31299

Capsulorhexis Forceps
short handle, extra delicate
cross serrated grasping tips
12 mm straight shafts, angled
overall length 83 mm





Geuder® 09 Forceps
Capsulorhexis

GS-32926

AUFFARTH
**MICS Double Cross Action
Capsulorhexis Forceps**
interlocking micro serrated tips
12 mm curved shafts, angled 130°
with overlap system
prevents inadvertent iris capture
for 1.6 mm incisions
overall length 100 mm









GS-32928


INAMURA
MICS Capsulorhexis Forceps
cross action
interlocking micro serrated tips
105 mm curved shafts tip to pivot
angled 130°
overall length 103 mm





GS-32935


H.-R. KOCH
MICS Capsule Forceps
short handle
22 Gauge/0.7 mm curved shaft with
micro grasping tips
overall length 98 mm





GS-32938

H.-R. KOCH
MICS Capsule Forceps
22 Gauge/0.7 mm curved shaft with
micro grasping tips
6 mm titanium handle, blue-colored
overall length 144 mm





Geuder® 09 Forceps
Micro Incision

32936

TOBIAS NEUHANN


Micro Forceps

multifunction asymmetrical tips with grasping platform

23 Gauge/ 0.6 mm straight shaft

6 mm titanium handle, blue-colored

overall length 155 mm



38244
38242

IKEDA

Capsulorhexis Forceps

interchangeable tips

23 Gauge/ 0.6 mm curved shaft

overall length 40 mm

textured, micrograsping tips "BEAK"

textured, micrograsping tips

for use with JAPAN MODEL Handle

38240



38240

JAPAN MODEL

Handle

for capsule scissors and forceps

interchangeable tips

6 mm, titanium

overall length 90 mm




32944

HATTENBACH

Capsulorhexis Forceps

25 Gauge/ 0.5 mm angled shaft with micro grasping tips

overall length 88 mm




38624

SZURMAN

DMEK Preparation Forceps

for removal of Descemet layer with smooth round tips to prevent tearing the donor endothelium

overall length 75 mm



Geuder® 09 Forceps
DMEK, DSAEK, LASIK

38628

YÖRÜK

DMEK Transplant Forceps

atraumatic grasping and removal of Descemet membrane

grasping area 6 mm

overall length 77 mm



38634

GERTEN

DSAEK/DMEK Rhexis Forceps

for Descemetorhexis

reversed triangular tips to score and remove the DM from the host cornea

23 Gauge/ 0.6 mm curved shaft

overall length 142 mm




33428

PFAEFFL

Flaporhexis Forceps

2.5 mm x 1.5 mm cross serrated tips, oval to lift Femto-LASIK flaps

overall length 103 mm




31899

LASIK Corneal Flap Forceps

angled shafts, 9 mm

highly polished round and smooth shafts

overall length 108 mm



18914

BURATTO

LASIK Corneal Flap Forceps

curved shafts

cross serrated tips

0.8 mm diameter

overall length 100 mm



Geuder® 09 Forceps
Lenticule

GS-33960

RUPAL SHAH

Lenticule Forceps


3.5 mm serrated tips

20 Gauge/ 0.9 mm curved shaft

to retrieve lenticule during SMILE procedures

6 mm titanium handle, blue-colored

overall length 144 mm



GS-33961

Lenticule Forceps


7.5 mm shafts

with fine grooves

angled to side

to retrieve lenticule during SMILE procedures

overall length 99 mm



Geuder® 09 Forceps
Nucleus Prechoppers, Implantation

GS-31518


Nucleus Spreading Forceps

cross-action

3.0 mm textured paddle tips

for cracking the nucleus

overall length 112 mm



GS-31519

AKAHOSHI

Nucleus Spreading Forceps

Phaco Prechopper

4.5 mm delicate sharp tips, 1 mm profile

2.0 mm cross-action

for cracking, dividing, and rotating the nucleus

overall length 105 mm



GS-31720

AKAHOSHI



Prechopper

for micro-incisions, straight shafts,

1.7 mm cross-action

round handle

overall length 125 mm



GS-33440

VON BARSEWISCH

Anterior Chamber Forceps

extra delicate

8.5 mm angled shafts

1.5 mm platform

overall length 72 mm




GS-32165

FECHNER

Implantation Forceps

smooth curved shafts

overall length 80 mm



Geuder® 09 Forceps
Nucleus Prechoppers, Implantation

GS-18935

DODEN
Implantation Forceps

curved
delicate lens holding shafts
overall length 110 mm



GS-31962

KOHNEN
Implantation Forceps

silicone tube coated
10 mm angled shafts
overall length 107 mm



Geuder® 09 Forceps
IOL

GS-18918

DICK
Fixation and Holding Forceps

for Artisan/Verisyse IOL
with guide mark for lens maximum hold
and stability
18 Gauge/ 12 mm curved shaft
overall length 143 mm



GS-18928

SEKUNDO
Implantation Forceps

for retropupillary implantation of
Iris-Claw-IOL
20 Gauge/ 0.9 mm curved shaft with
lens holding plate, three-point fixation
and IOL centration guide mark
6 mm titanium handle, blue-colored
overall length 143 mm



Geuder® 09 Forceps
IOL

GS-31902

Folder

foldable IOL insertion forceps
cross action
65 mm concave jaws
overall length 103 mm



GS-31903

SPALECK
Folder

foldable IOL insertion forceps
for silicone and acrylic IOLs
angled shafts, concave jaws
overall length 102 mm



GS-31908

Folder

IOL folding forceps
for easy pickup and folding acrylic IOLs
overall length 105 mm



GS-31994

H.-R. KOCH
Prefolder

for foldable silicone, hydrogel or
acrylic IOLs
overall length 135 mm



GS-31905

H.-R. KOCH
Folder


for foldable 3-piece IOLs
6.0 mm concave jaws
overall length 92 mm



Geuder® 09 Forceps
IOL

GS-31901


Holding Forceps
for foldable IOLs
highly polished round and smooth shafts
overall length 108 mm





GS-18941

STROBEL Holding Forceps
for Pharmacia CeeOn Edge IOL
overall length 110 mm





GS-32823


Explantation Forceps
for soft IOLs
2.0 mm x 1.4 mm paddle tip with
interlocking serrations and guide pin
overall length 108 mm





GS-32997

SNYDER - OSHER IOL Forceps
to hold the IOL while cutting it
inside the capsular bag
2 mm long jaws with
0.5 mm serrated tips
22 Gauge/ 0.7 mm lightly curved shaft
overall length 99 mm






Geuder® 09 Forceps
Superior Rectus

GS-18820


REMKY Superior Rectus Forceps
angled
1 x 2 teeth
overall length 109 mm





GS-18800


ELSCHNIG Superior Rectus Forceps
standard, straight
1 x 2 teeth
overall length 110 mm





GS-19105


LEYDHECKER Superior Rectus Forceps
curved, tip width 0.8 mm
1 x 2 teeth, 1.1 mm
overall length 110 mm





GS-32335

CASTROVIEJO Superior Rectus Forceps
angled with tying platform
1 x 2 teeth, 0.9 mm
overall length 100 mm





Geuder® 09 Forceps
Muscle, Nasal, Cerclage, Retinal

GS-19112

BANGERTER
Muscle Forceps
10 teeth
self retaining
overall length 100 mm



GS-19120
GS-19125

SATTLER
Muscle Forceps
with slide lock
overall length 100 mm
3 teeth, bulbous tip, right
3 teeth, bulbous tip, left



GS-32350

TROELTSCH
Nasal Forceps
serrated tips
for packing and dressing during ENT procedures
overall length 150 mm



GS-19107

ESSEN MODEL
Cerclage Forceps
strongly curved
3 mm platform with dual grooves
overall length 104 mm



GS-33434


PFENNIGWERTH
Holding Forceps
for vitreoretinal lenses
curved with concave lens holding tips
overall length 110 mm



Geuder® 09 Forceps
Muscle, Nasal, Cerclage, Retinal

GS-33438


Forceps
for HEIMANN retinal tacks
overall length 100 mm



Geuder® 09 Forceps
Scleral Plug, Trocar

GS-33439

Scleral Plug Forceps
cross action handle, angled jaws
to securely hold plugs during insertion
and removal
overall length 105 mm



GS-33752

WATZKE
Spreading Forceps
to stretch silicone sleeves
cross action handle and
serrated tips
overall length 115 mm




GS-37430

Trocar Forceps
for trocar systems with groove
overall length 110 mm






G-33831
REIMANN
Vitreoretinal Forceps
for foreign body removal
with retractable magnet and 3-prongs
14 Gauge/ 2.1 mm straight shaft
shaft length 35 mm
overall length 138 mm




Vitreoretinal Forceps
6 mm titanium handle
blue-colored
overall length 155 mm


G-36202
straight
serrated tips
20 Gauge/ 0.9 mm




G-36212
straight
silicone coated
20 Gauge/ 0.9 mm





G-36256
G-36342
G-36322
straight
serrated "crocodile"
tips
20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm





G-36262
foreign body
forceps
20 Gauge/ 0.9 mm





G-36362
G-36366
TANO
straight, serrated tips
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm



G-36216
G-36286
straight
grip tips
20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm



G-36292
G-36282
G-36312
straight
endgripping,
extra delicate
0.2 mm x 0.05 mm
20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm



Single-Use Interchangeable Forceps
available on page [176](#)



Geuder® 09 Forceps
Chalazion, Ptosis

⌘-18680

LAMBERT
Chalazion Forceps
round, 9 mm inside diameter
solid and fenestrated plates
overall length 92 mm

⌘-18685

⌘-18690

DESMARRES
Chalazion Forceps
oval
solid and fenestrated plates
20 mm, overall length 88 mm
25 mm, overall length 90 mm

⌘-18700

CAUER
Chalazion Forceps
oval, 13 mm x 11 mm inside diameter
interlacing serrated tips
solid and fenestrated plates
overall length 99 mm

⌘-18715

BERKE
Ptosis Forceps
22 mm long serrated tips
slide lock mechanism
overall length 110 mm

Geuder® 09 Forceps
Mosquito

⌘-18170

⌘-18171

HARTMANN
Mosquito Forceps
serrated
18 mm jaws
overall length 100 mm
straight
curved

⌘-18180

⌘-18181

HALSTEAD
Mosquito Forceps
serrated
25 mm jaws
overall length 125 mm
straight
curved

⌘-18185

⌘-18182

HALSTEAD
Mosquito Forceps
1 x 2 teeth
serrated
25 mm jaws
overall length 125 mm
straight
curved

Geuder® 09 Forceps
Towel, Serrefine

⌘-18160	BACKHAUS Towel Forceps overall length 90 mm	
⌘-33200	SCHAEDEL Towel Forceps overall length 90 mm	
⌘-18190	DIEFFENBACH Serrefine straight, overall length 35 mm	
⌘-18195	DIEFFENBACH Serrefine curved, overall length 35 mm	
⌘-31420	DIEFFENBACH Serrefine straight, overall length 50 mm	
⌘-18197	DIEFFENBACH Serrefine extra delicate, overall length 50 mm	

Geuder® 10 Forceps