

Geuder 10
Needle Holders

Nadel- und
Klingenhalter



Geuder 10 Needle Holders

G-17500
G-17501

BARRAQUER
Needle Holder
10 mm delicate jaws
0.8 x 0.55 mm tip
overall length 120 mm
curved, without lock
curved, with lock




G-17505
G-17506

BARRAQUER
Needle Holder
10 mm delicate jaws
0.8 x 0.45 mm tip
overall length 120 mm
curved, without lock
curved, with lock



G-17450
G-17452

JACOBI
Needle Holder
5 mm extra delicate jaws
1.0 x 0.5 mm tips
overall length 115 mm
curved, without lock
curved, with lock



G-17460
G-17462

BARRAQUER
Needle Holder
7.5 mm extra delicate jaws
0.6 x 0.3 mm tip
overall length 105 mm
curved, without lock
curved, with lock



G-17535
G-17530

BARRAQUER - TROUTMAN
Needle Holder
10 mm extra delicate jaws
0.8 x 0.6 mm tip
handle diameter 6 mm
overall length 100 mm
curved, without lock
curved, with lock





BARRAQUER
Needle Holder
11 mm delicate jaws
0.9 x 0.8 mm tips
handle diameter 8 mm
overall length 140 mm
curved, without lock

G-17510
G-17511



curved, with lock



BARRAQUER
Needle Holder
11 mm extra delicate jaws
0.9 x 0.55 mm tips
handle diameter 8 mm
overall length 140 mm
curved, without lock

G-17515
G-17516



curved, with lock



BARRAQUER
Needle Holder
11 mm delicate jaws
1.1 x 0.75 mm tips
handle diameter 10 mm
curved, without lock
overall length 140 mm

G-17520



curved, with lock



G-32315

BARRAQUER
Needle Holder
curved, without lock
8 mm standard jaws
1.2 x 0.8 mm tips
original model with longitudinal groove
handle diameter 10 mm
overall length 130 mm



curved, with lock



BARRAQUER
Needle Holder
12 mm standard jaws
1.6 x 1.0 mm tips
handle diameter 10 mm
overall length 140 mm
curved, without lock

G-32320
G-32321



curved, with lock



CASTROVIEJO
Needle Holder
10.5 mm delicate jaws
1.5 x 0.8 mm tips
overall length 140 mm
straight, without lock

G-17560
G-17565



straight, with lock




CASTROVIEJO
Needle Holder
10.5 mm delicate jaws
1.5 x 0.8 mm tips
overall length 140 mm
curved, without lock
curved, with lock

17570
17575



Combined Needle Holder/ Scissors
curved 12 mm blades with
3 mm needle holder and tying platform
used to hold, place and cut fine suture material
overall length 125 mm



17440



HATTENBACH
Hybrid Needle Holder
angled, for the anterior chamber
25 Gauge/ 0.5 mm
overall length 90 mm

for use in combination with
HATTENBACH Hybrid Tying Forceps 32947



32949



KALT
Needle Holder
standard jaws
serrated, with longitudinal groove
thumb rest and catch
overall length 135 mm
straight, cross serrated jaws
curved, cross serrated jaws

17610
17615



ARRUGA
Needle Holder
straight
one jaw smooth other serrated
14 mm heavy jaws
thumb rest and catch
overall length 140 mm

17640



HEIDELBERG MODEL
Needle Holder
straight
14 mm heavy jaws
with ratchet lock
overall length 125 mm

17660

☒-19320

CASTROVIEJO

Razor Blade Holder

6 mm round handle

overall length 120 mm



☒-19330

TROUTMAN

Razor Blade Holder

straight

6 mm round handle

overall length 95 mm



☒-18150

Holder for Absorbent Strips/ Swabs

serrated

to secure Sterile Absorbent Swabs ☒-18127

overall length 175 mm