

# Scleral Depressors, Localizers

Lokalisatoren,  
Skleradepressoren



GS-32152

LAQUA

Scleral Depressor

with pocket clip

3 x 1.8 mm teardrop-shaped tip

overall length 128 mm



GS-32150

LAQUA


Scleral Depressor

double-ended with pocket clip

3 x 1.8 mm teardrop-shaped tip

2 mm diameter scleral marker

overall length 134 mm



GS-32153

LUCKE

Scleral Depressor

with pocket clip

3.8 x 2.0 mm teardrop-shaped tip

overall length 128 mm



GS-32027

KREISSIG


Scleral Depressor

double-ended with pocket clip

conical tip with 2.4 mm diameter

and 6.0 mm bar

overall length 128 mm



GS-15710

MEYER - SCHWICKERATH

Localizer and Scleral Depressor

for fixation and rotation



GS-32715

SCHOCKET


Scleral Depressor

double-ended with pocket clip

3 x 2 mm teardrop tip and

6.0 mm bar

overall length 132 mm



GS-32714

Scleral Depressor

double-ended

2 mm and 4 mm diameter

diamond dusted tips

overall length 120 mm



GS-37445

HATTENBACH

Triangle Scleral Depressor

for gentle buckling and firm grip without the risk of rupture

overall length 133 mm

