

Geuder[®]

13

Probes,
Dilators

Sonden,
Dilatatoren

Lacrimal Probes 149

Pigtail, Trabecular, U-Shaped 151



Geuder[®]

13 Probes, Dilators

Lacrimal Probes

G-150301	WILDER Lacrimal Probe delicate, long taper 0.45 mm tip diameter overall length 104 mm	
G-150302	WILDER Lacrimal Probe medium, medium taper 0.50 mm tip diameter overall length 104 mm	
G-150303	WILDER Lacrimal Probe strong, short taper 0.55 mm tip diameter overall length 104 mm	
G-150401	WILDER Lacrimal Probe slim delicate, long taper 0.30 mm tip diameter overall length 95 mm	
G-150402	WILDER Lacrimal Probe slim medium, medium taper 0.35 mm tip diameter overall length 95 mm	
G-150403	WILDER Lacrimal Probe slim strong, short taper 0.40 mm tip diameter overall length 95 mm	



Geuder® 13 Probes, Dilators

Lacrimal Probes

GS-15020	GIESSEN MODEL Lacrimal Probe 11 mm tapered shaft, angled 0.4 mm tip width overall length 77 mm	
GS-150101	BOWMAN Lacrimal Probe cylindrical overall length 130 mm size 0000/000	
GS-150102	size 00/0	
GS-150103	size 1/2	
GS-150104	size 3/4	
GS-150105	size 5/6	
GS-150151	BOWMAN Lacrimal Probe bulbous tips overall length 130 mm size 0000/000	
GS-150152	size 00/0	
GS-150153	size 1/2	
GS-150154	size 3/4	
GS-150155	size 5/6	
GS-33193	JUENEMANN Lacrimal Probe with 72 mm probe and 18 mm grip plate to insert suture material provided with Cleaning Adapter GS-31022 0.6 mm diameter	 
GS-33194	0.8 mm diameter	
GS-33195	1.0 mm diameter	
GS-15050	KELLNAR Retrograde Probe for canalicular repair 0.5 mm suture hole overall length 95 mm right	
GS-15055	left	

Geuder® 13 Probes, Dilators

Pigtail, Trabecular, U-Shaped

GS-15095	Pigtail Probe double-ended for canalicular repair with 0.4 mm suture holes overall length 140 mm	  
GS-15096	Pigtail Probe double-ended for canalicular repair with 1.0 mm suture holes overall length 140 mm	  
GS-15125	HARMS Trabecular Probe dual 9.5 mm curved probes 3 mm spread overall length 46 mm right	
GS-15130	left	
GS-15136	THOMAS NEUHANN U-Shaped Probe dual 9 mm curved probes 3 mm spread overall length 12 mm 0.30 mm probe diameter	
GS-15138	0.25 mm probe diameter	
GS-15137	THOMAS NEUHANN Holder with lock for U-shaped Probes GS-15136 and GS-15138 overall length 92 mm	