

Geuder 14
Cannulas,
I/A Instruments







Kanülen,
Spül- und
Sauginstrumente (I/A)






- Lacrimal 153
- Cyclodialysis 155
- Hydrodissection 155
- Irrigation Vectis 156
- Viscoelastic 156
- Capsule Polishers 158
- LASIK 159
- Bimanual Irrigation/ Aspiration 160
- DMEK/ DSAEK Irrigation 162
- Retinal Backflush 165
- Retinal Aspiration, PFCL, etc. 167
- Miscellaneous 169

Geuder 14 Cannulas, I/A Instruments
Lacrimal

G-33135	Lacrimal Cannula conical, straight tapers from 19 to 25 Gauge 25 mm shaft stainless steel		
G-15141	Lacrimal Cannula conical, strongly curved tapers from 19 to 25 Gauge 25 mm shaft stainless steel		
G-15151	Lacrimal Cannula conical, angled at 10 mm tapers from 19 to 25 Gauge 25 mm shaft stainless steel		
G-15161	Lacrimal Cannula conical, angled at 18 mm tapers from 19 to 25 Gauge 25 mm shaft stainless steel		
G-15190	ANEL Lacrimal Cannula straight 15 mm shaft 23 Gauge/ 0.64 mm		
G-15195	ANEL Lacrimal Cannula strongly curved 16 mm shaft 23 Gauge/ 0.64 mm		
G-33160	Lacrimal Cannula cylindrical, straight 20 Gauge/ 0.9 x 35 mm		



GS-15170	BONN MODEL Lacrimonal Cannula with probe tip 1.15 mm diameter olive tip		
GS-15271 GS-15272 GS-15273	BANGERTER Lacrimonal Probe Cannula straight, 0.6 mm side port 23 Gauge/ 0.64 x 60 mm 21 Gauge/ 0.80 x 40 mm 23 Gauge/ 0.64 x 40 mm		
GS-33197	JUENEMANN Lacrimonal Probe Cannula to insert suture material 21 Gauge/ 0.8 x 35 mm		
GS-33198	JUENEMANN Lacrimonal Probe Cannula with grip plate 65 mm probe to insert suture material 21 Gauge/ 0.8 x 65 mm		
GS-15377	Retrobulbar Injection Cannula 25 Gauge/ 0.5 x 40 mm		
GS-15200	CRANZ Eye Douche irrigation through five front ports silver		

GS-15275	Cyclodialysis Spatula Cannula 12 mm long spatula anterior 0.2 mm port 26 Gauge/ 0.45 mm		
GS-31355	BUSSE Cyclodialysis Spatula Cannula 12 mm long spatula dual 0.2 mm port 26 Gauge/ 0.45 mm		

GS-15250 GS-32166 GS-32167	SAUTTER Hydrodissection Cannula front opening flattened tip vertical, 27 Gauge/ 0.4 mm vertical, 26 Gauge/ 0.45 mm horizontal, 27 Gauge/ 0.4 mm	 vertical GS-15250, GS-32166	 horizontal GS-32167	
GS-32155	TUEBINGEN MODEL Hydrodissection Cannula 30 Gauge/ 0.3 mm			
GS-15185	FRANKFURT MODEL Hydrodissection Cannula 27 Gauge/ 0.4 mm			
GS-31895	Hydrodissection Cannula vertically flattened tip front opening 24 Gauge/ 0.6 x 20 mm			




Geuder® 14 Cannulas, I/A Instruments

Irrigation Vectis

GS-16115	Irrigation Vectis Cannula 4 front ports 5 mm wide lens loop curved		
GS-32759	Irrigation Vectis Cannula small model 4 front ports 4 mm wide lens loop curved		




Geuder® 14 Cannulas, I/A Instruments








Viscoelastic

GS-31807	Healon - Viscoelastic Cannula flattened tip, angled at 6 mm 30 mm shaft 21 Gauge/ 0.8 mm		
GS-22619	PH. JACOBI Trabecular Aspiration Cannula angled at 13 mm 24 Gauge/ 0.6 mm		
GS-15188 GS-15189	BURK Irrigation Trabeculotomy Cannula 8 mm curved irrigation probe 30 Gauge/ 0.3 mm angled right angled left		
GS-33199	NAGATA Sub-Tenon Injection Cannula for atraumatic penetration and injection of anesthetic agent 24 Gauge/ 0.6 mm		



Geuder® 14 Cannulas, I/A Instruments

Irrigation/ Aspiration

GS-15180	Irrigation Cannula bulbous tip, 0.75 mm diameter tip extra delicate 26 Gauge/ 0.45 mm		
GS-15230 GS-15240	BONN MODEL Irrigation Cannula flattened tip, front opening delicate, 0.5 x 1.5 mm medium, 0.6 x 1.9 mm		
GS-32760 GS-32765	MUELLER Irrigation Cannula J-shaped with probe tip 26 Gauge/ 0.45 x 19 mm angled left angled right		
GS-16138	RENTSCH Irrigation Cannula J-shaped angled at 10 mm 26 Gauge/ 0.45 x 24 mm		
GS-16106	SIMCOE Irrigation/ Aspiration Cannula dual shaft 0.3 mm diameter aspiration port front opening irrigation port 23 Gauge/ 0.64 mm with silicone tube and Luer-Lock hub		
GS-15225 GS-15226	Aspiration Cannula lightly curved, front opening delicate, 19 Gauge/ 1.1 mm standard, 17 Gauge/ 1.5 mm		

GS-16125	KRATZ Capsule Polishing Cannula angled at 2 mm half-textured tip 23 Gauge/ 0.64 mm		
GS-16133	TEPING Aspiration Polishing Cannula 0.3 mm posterior port half-textured tip, curved 23 Gauge/ 0.64 mm		
GS-16130	Capsule Polishing Cannula with half-textured ball tip 18 Gauge/ 1.2 mm		
GS-16132	PAPST Capsule Polishing Cannula 0.75 mm hemisphere textured tip 25 Gauge/ 0.5 mm		
GS-32529	KYU Capsule Polishing Cannula 2 mm disc-shaped tip textured 22 Gauge/ 0.7 mm		
GS-32528	GREITE Capsule Polishing Cannula curette tip textured 22 Gauge/ 0.7 mm		
GS-32774	BRAUWEILER Aspiration Cannula angled at 13 mm textured polisher tip 0.3 mm beveled port 21 Gauge/ 0.8 mm		

GS-32769	H.-R. KOCH Aspiration Cannula angled at 14 mm 0.3 mm anterior port half-textured, spatulated tip 21 Gauge/ 0.8 mm		
GS-32776	BRAUWEILER Irrigation Cannula angled at 12 mm front opening 21 Gauge/ 0.8 mm		
GS-32772	BRAUWEILER Irrigation Cannula angled at 13 mm textured 0.8 mm iris retractor and capsule polisher 21 Gauge/ 0.8 mm		

GS-33909	W. BAUER LASIK Cannula double-armed 6 ports for irrigation of the flap and stromal bed 25 Gauge/ 0.5 mm		
GS-33911	ZIRM LASIK Irrigation Cannula flattened and textured tip with 5 irrigation ports 20 Gauge/ 0.9 mm		

GS-15310

KRUMEICH


Corneal Irrigation Cannula

for big bubble technique during DALK

semi-sharp dissector tip

0.2 mm posterior port

27 Gauge/ 0.4 mm





GS-32015

Irrigation Cannula Handpiece

male/ female connectors

overall length 112 mm



GS-32771

BRAUWEILER

Irrigation Cannula Handpiece

male/ female connectors

overall length 78 mm



GS-32773

BRAUWEILER

Irrigation Cannula Handpiece

male/ male connectors

overall length 80 mm



GS-22632

SCHMACK

Irrigation Handpiece

21 Gauge/ 0.8 mm

0.5 mm port



GS-22633

SCHMACK

Aspiration Handpiece

21 Gauge/ 0.8 mm

0.3 mm port



GS-22100


D.D. KOCH - KOHNEN

Aspiration Handpiece

0.3 mm aspiration port

textured tip

22 Gauge/ 0.7 mm





GS-22102

23 Gauge/ 0.6 mm

GS-22101

D.D. KOCH - KOHNEN

Irrigation Handpiece

front opening

22 Gauge/ 0.7 mm






GS-22106

Aspiration Handpiece

0.25 mm aspiration port

textured, conical tip

oval tube 0.95 x 0.55 mm






GS-22107


Irrigation Handpiece

0.5 mm irrigation port

conical tip

oval tube 0.95 x 0.55 mm








GS-32011

Aspiration Handpiece

0.3 mm aspiration port

textured, conical tip

21 Gauge/ 0.8 mm





GS-32012

Irrigation Handpiece

dual 0.5 mm irrigation ports


textured, conical tip

21 Gauge/ 0.8 mm












Geuder® 14 Cannulas, I/A Instruments
Bimanual Irrigation/ Aspiration




GS-32018	BRAUWEILER Aspiration Handpiece angled at 13 mm 0.3 mm beveled port textured tip 21 Gauge/ 0.8 mm	 
GS-32019	H.-R. KOCH Aspiration Handpiece angled at 14 mm 0.3 mm anterior port half-textured, spatulated tip 21 Gauge/ 0.8 mm	 
GS-32016	BRAUWEILER Irrigation Handpiece angled at 12 mm front opening 21 Gauge/ 0.8 mm	 


Geuder® 14 Cannulas, I/A Instruments
DMEK/ DSAEK Irrigation

GS-38602	THOMAS NEUHANN Irrigation Descemet Scraper scraper width 0.8 mm scraper length 0.6 mm rhombic 0.9 x 0.6 mm irrigation tube for peripheral scraping of the recipient bed overall length 94 mm	 
GS-38632	SZURMAN DMEK Irrigation Handpiece curved, front opening tube 23 Gauge/ 0.6 mm	 

Geuder® 14 Cannulas, I/A Instruments
Coaxial Irrigation/ Aspiration








		  
GS-24318 GS-24418	PEP 2.2 mm I/A Handpiece Pure Efficiency Phaco for 2.2 mm incisions 0.3 mm aspiration port coaxial I/A straight tip coaxial I/A 30° angled tip with GS-24311 PEP Infusion Sleeve	straight GS-24318 angled GS-24418
GS-24311	PEP 2.2 mm Infusion Sleeve Pure Efficiency Phaco for 2.2 mm incisions purple translucent silicone	 
		  
GS-24332 GS-24432	PEP 2.4 mm I/A Handpiece Pure Efficiency Phaco for 2.4 mm incisions 0.3 mm aspiration port coaxial I/A straight tip coaxial I/A 30° angled tip with GS-24314 PEP Infusion Sleeve	straight GS-24332 angled GS-24432
GS-24314	PEP 2.4 mm Infusion Sleeve Pure Efficiency Phaco for 2.4 mm incisions green translucent silicone	 

		
I/A Handpiece coaxial 30° angled tip		
┐-24430	0.3 mm aspiration port	
┐-24435	0.35 mm aspiration port	
┐-24440	0.4 mm aspiration port with ┐-24100 Infusion Sleeve	
┐-24100	PEP 2.4 mm Infusion Sleeve for I/A handpieces blue silicone	
		
PAPST I/A Handpiece coaxial		
infusion tube 16 Gauge/ 1.6 mm		
┐-22525	aspiration port 0.3 mm, straight	
	aspiration port 0.3 mm, angled 30°	
┐-22540		

┐-33045	Backflush Handpiece with Luer cone overall length 110 mm	
┐-31805	Backflush Handpiece with silicone chamber and Luer cone overall length 110 mm	
┐-37002	Backflush Handpiece with silicone chamber and Luer-Lock connector overall length 115 mm	
┐-37000	Backflush Handpiece with silicone chamber and chuck overall length 110 mm	
for backflush tips ┐-37010 up to ┐-37080		
Silicone Chamber for backflush handpieces ┐-37000, ┐-37002, ┐-31805		
┐-31806	reusable	
┐-34806	single use, 10 pcs. per box, sterile	
┐-37005	Backflush Handpiece with large silicone chamber rotatable chuck	
for backflush cannulas ┐-37010 up to ┐-37080		
┐-37006	Silicone Chamber large, for Backflush Handpiece ┐-37005	
Backflush Cannula 35 mm shaft		
┐-33050	20 Gauge/ 0.9 mm	
┐-33054	23 Gauge/ 0.6 mm	
┐-33053	25 Gauge/ 0.5 mm	
		







Geuder® 14 Cannulas, I/A Instruments

Retinal Backflush

GS-33070	Backflush Cannula 35 mm shaft 20 Gauge/ 0.9 mm shaft 25 Gauge/ 0.5 mm tip 3.5 mm tapered tip		
GS-33073	KAMPIK Backflush Cannula 35 mm shaft inside diameter 22 Gauge/ 0.7 mm outside diameter 20 Gauge/ 0.9 mm		
GS-33075	Backflush Cannula 35 mm shaft 0.8 mm silicone brush 20 Gauge/ 0.9 mm		
GS-37010 GS-37011 GS-37012	Backflush Tip 35 mm shaft 20 Gauge/ 0.9 mm 23 Gauge/ 0.6 mm 25 Gauge/ 0.5 mm		
GS-37020	Backflush Tip 35 mm shaft 0.8 mm silicone brush 20 Gauge/ 0.9 mm		
GS-37035	Backflush Tip 35 mm shaft 0.8 mm silicone tip 20 Gauge/ 0.9 mm		
GS-37080	Backflush Tip 35 mm shaft 20 Gauge/ 0.9 mm shaft 23 Gauge/ 0.6 mm tip 2 mm tapered tip overall length 35 mm		

Geuder® 14 Cannulas, I/A Instruments

Retinal Aspiration, PFCL, etc.


GS-33051 GS-33052	WEBER Aspiration Cannula for viscous fluids inside diameter 22 Gauge/ 0.7 mm outside diameter 20 Gauge/ 0.9 mm 5 mm shaft 10 mm shaft		
GS-37390	Aspiration Cannula 23 Gauge/ 0.6 mm 10 mm shaft		
GS-33056	ROIDER Aspiration Cannula for viscous fluids 0.7 mm side port 19 Gauge/ 1.0 mm		
GS-33055	PFCL Cannula with dual shaft 24 Gauge/ 0.55 mm 25 Gauge/ 0.5 mm, lateral tube		
GS-33057	CHANG PFCL Cannula for injection of heavy fluids dual bore, coaxial 25 Gauge/ 0.5 mm tip 20 Gauge/ 0.9 mm shaft		
GS-32699 GS-32698	HEIDELBERG MODEL Cannula for injection or aspiration of viscous fluids and Densiron® 68, bevel 30° 19 Gauge/ 1.1 mm 18 Gauge/ 1.2 mm		
GS-33060 GS-33065 GS-33066	Stopper for viscous fluid aspiration with tube connection for single-use syringe for 5 ml syringe for 10 ml syringe for 20 ml syringe		

GS-32697 **Pressure Tube**
for injection of viscous fluids
Luer-Lock female/ male
reusable



GS-30745 **Pressure Connection Tube**
for megaTRON S3/ S4 HPS
for Holder GS-28764 for pneumatic silicone oil injection
length 2 m

GS-28764 **Adapter**
with protective cover for glass syringes
to inject viscous fluid pneumatically



GS-19961 **Syringe Clip**
customized labeling
for 5 ml syringes



GS-32665 **Silicone Bulb**
with Luer cone
20 ml
GS-32670 60 ml

