

Geuder®

15

Single-Use Instruments

Einmalinstrumente

- Trephines and Punches

171
- Surgical Knives

173
- Scissors and Forceps

176
- Laser Probes

178
- Fiber Optics

180
- Vitreous Cutters

181
- Trocar Systems

181
- DMEK

182
- Cannulas

183
- Injection Syringes

187
- Filter, Stop Cock, Injection Tube

188
- Oil Injection

189
- Cautery, Scalpel Handle, Absorbent Swabs, Marking Pen

190
- Silicone Tires

191
- Silicone Circling Bands

192
- Silicone Strips

192
- Silicone Sleeves

193
- Silicone Sponge Implants

194



Geuder®

15 Single-Use Instruments

Trephines and Punches

NanoEdge

NanoEdge Trephines
length 16 mm
sterile, single-use

GS-520400	4.00 mm	GS-520825	8.25 mm	GS-521300	13.00 mm
GS-520550	5.50 mm	GS-520850	8.50 mm	GS-521350	13.50 mm
GS-520600	6.00 mm	GS-520875	8.75 mm	GS-521400	14.00 mm
GS-520625	6.25 mm	GS-520900	9.00 mm	GS-521450	14.50 mm
GS-520650	6.50 mm	GS-520950	9.50 mm	GS-521500	15.00 mm
GS-520675	6.75 mm	GS-521000	10.00 mm	GS-521550	15.50 mm
GS-520700	7.00 mm	GS-521050	10.50 mm	GS-521600	16.00 mm
GS-520725	7.25 mm	GS-521100	11.00 mm	GS-521650	16.50 mm
GS-520750	7.50 mm	GS-521150	11.50 mm	GS-521700	17.00 mm
GS-520775	7.75 mm	GS-521200	12.00 mm	GS-521800	18.00 mm
GS-520800	8.00 mm	GS-521250	12.50 mm		

NanoEdge



- GS-500001

NanoEdge Single-Use Corneal Cutting Block
Teflon, sterile
white-colored
outside diameter 18.8 mm
height 9.5 mm

NanoEdge




NanoEdge Trephine Punches
sterile, single-use

GS-540600	6.00 mm	GS-540800	8.00 mm
GS-540625	6.25 mm	GS-540825	8.25 mm
GS-540650	6.50 mm	GS-540850	8.50 mm
GS-540675	6.75 mm	GS-540875	8.75 mm
GS-540700	7.00 mm	GS-540900	9.00 mm
GS-540725	7.25 mm	GS-540925	9.25 mm
GS-540750	7.50 mm	GS-540950	9.50 mm
GS-540775	7.75 mm		

NanoEdge






NanoEdge Vacuum Trephine Punches
sterile, single-use

GS-560600	6.00 mm	GS-560725	7.25 mm	GS-560850	8.50 mm
GS-560625	6.25 mm	GS-560750	7.50 mm	GS-560875	8.75 mm
GS-560650	6.50 mm	GS-560775	7.75 mm	GS-560900	9.00 mm
GS-560675	6.75 mm	GS-560800	8.00 mm	GS-560925	9.25 mm
GS-560700	7.00 mm	GS-560825	8.25 mm	GS-560950	9.50 mm



NanoEdge




NanoEdge Vacuum Trephine radial
sterile, single-use

GS-570600	6.00 mm	GS-570750	7.50 mm	GS-570850	8.50 mm
GS-570650	6.50 mm	GS-570775	7.75 mm	GS-570875	8.75 mm
GS-570700	7.00 mm	GS-570800	8.00 mm	GS-570900	9.00 mm
GS-570725	7.25 mm	GS-570825	8.25 mm		



NanoEdge




protect, with sliding guard

NanoEdge Paracentesis Knife

straight
6 pcs. per box

GS-34026
GS-34064

15°
15° protect



GS-34031

22,5°



NanoEdge


GS-34036
GS-34065

30°
30° protect



GS-34041

45°



NanoEdge Stab Blade Tip

straight
6 pcs. per box

GS-34042
GS-34043
GS-34044

15°
30°
45°



NanoEdge




NanoEdge MVR Stiletto

straight
6 pcs. per box

GS-34176
GS-34181
GS-34192
GS-34198



diameter insertable
19 Gauge/ 1.1 mm
20 Gauge/ 0.9 mm
23 Gauge/ 0.6 mm
25 Gauge/ 0.5 mm

incision/ blade width
1.4 mm
1.2 mm
0.8 mm
0.7 mm



NanoEdge



G-42361 G-42561	Single-Use Vitreous Scissors UNO Colorline straight, sterile			
	23 G, 1 pc. per box 25 G, 1 pc. per box			
G-46362 G-46562	Single-Use Vitreous Scissors UNO Colorline curved, sterile			
	23 G, 6 pcs. per box 25 G, 6 pcs. per box			
G-46341 G-46541	Single-Use Vitreous Forceps Endgripping UNO Colorline straight, extra delicate, endgripping, sterile			
	23 G, 6 pcs. per box 25 G, 6 pcs. per box			
G-46342 G-46542	HATTENBACH Single-Use Vitreous Forceps Wave Design UNO Colorline straight, wave design, sterile			
	0.2 mm x 0.5 mm 23 G, 6 pcs. per box 25 G, 6 pcs. per box			
G-42043 G-46343	Single-Use Vitreous Forceps Crocodile Design UNO Colorline straight, crocodile design, sterile			
	0.3 mm x 2.5 mm 20 G, 1 pc. per box 0.2 mm x 1.5 mm 23 G, 6 pcs. per box			
G-42544	Single-Use Tano Forceps UNO Colorline straight, Tano design, sterile			
	0.5 mm 25 G, 1 pc. per box			


G -40006 G -40008	Handle UNO Colorline 20 G, 23 G, 25 G titanium handle for UNO Colorline single-use tips 6 mm diameter 8 mm diameter				
G -34502	AMON Single-Use Forceps Needle for intrascleral IOL-fixation 27 Gauge/ 0.4 mm 6 pcs. per box, sterile for use with JAPAN MODEL Handle G -38240				 
G -38240	JAPAN MODEL Handle for capsule scissors and forceps interchangeable tips 6 mm, titanium overall length 90 mm				

CE pending



Straight Laser Probes UNO Colorline
Laser Connector: Geuder, Zeiss, Alcon

- 61508 20 G, 10 pcs. per box
- 61538 23 G, 10 pcs. per box
- 61558 25 G, 10 pcs. per box
- 61578 27 G, 10 pcs. per box




Straight Laser Probes UNO Colorline
Laser Connector: Iridex®

- 61601 20 G, 10 pcs. per box
- 61631 23 G, 10 pcs. per box
- 61651 25 G, 10 pcs. per box
- 61671 27 G, 10 pcs. per box




Curved Laser Probes UNO Colorline
Laser Connector: Geuder, Zeiss, Alcon

- 61509 20 G, 10 pcs. per box
- 61539 23 G, 10 pcs. per box
- 61559 25 G, 10 pcs. per box
- 61579 27 G, 10 pcs. per box





Curved Laser Probes UNO Colorline
Laser Connector: Iridex®

- 61603 20 G, 10 pcs. per box
- 61633 23 G, 10 pcs. per box
- 61653 25 G, 10 pcs. per box
- 61673 27 G, 10 pcs. per box




Bending Laser Probes UNO Colorline
Laser Connector: Geuder, Zeiss, Alcon

- 61501 20 G, 10 pcs. per box
- 61531 23 G, 10 pcs. per box
- 61551 25 G, 10 pcs. per box
- 61571 27 G, 10 pcs. per box



Bending Laser Probes UNO Colorline
Laser Connector: Index®

- 61503 20 G, 10 pcs. per box
- 61533 23 G, 10 pcs. per box
- 61553 25 G, 10 pcs. per box
- 61573 27 G, 10 pcs. per box




Illuminated Bending Laser Probes UNO Colorline
Light Connector: Geuder
Laser Connector: Geuder, Zeiss, Alcon

- 61505 20 G, 10 pcs. per box
- 61535 23 G, 10 pcs. per box
- 61555 25 G, 10 pcs. per box
- 61575 27 G, 10 pcs. per box



light connector




Illuminated Bending Laser Probes UNO Colorline
Light Connector: Alcon
Laser Connector: Geuder, Zeiss, Alcon

- 61506 20 G, 10 pcs. per box
- 61536 23 G, 10 pcs. per box
- 61556 25 G, 10 pcs. per box
- 61576 27 G, 10 pcs. per box




light connector



Geuder® 15 Single-Use Instruments
Fiber Optics

Single-Use Fiber Optic S
UNO Colorline
spot, sterile

☒ -46021 20 G, 6 pcs. per box
☒ -46321 23 G, 6 pcs. per box
☒ -46521 25 G, 6 pcs. per box



Single-Use Fiber Optic W
UNO Colorline
wide angle, sterile

☒ -46022 20 G, 6 pcs. per box
☒ -46322 23 G, 6 pcs. per box
☒ -46522 25 G, 6 pcs. per box



Single-Use Fiber Optic WS
UNO Colorline
wide angle, shielded, sterile

☒ -46023 20 G, 6 pcs. per box
☒ -46323 23 G, 6 pcs. per box
☒ -46523 25 G, 6 pcs. per box



HATTENBACH
Single-Use Flat Head Trocar
UNO Colorline
1 port for one-step technique

☒ -S03503 25 G, 6 pcs. per box



Geuder® 15 Single-Use Instruments
Vitreous Cutters



Single-Use High Speed Vitreous Cutter
UNO Colorline
single-blade, pneumatic,
with Geuder safety-connector, sterile

☒ -42001 20 G, 1 pc. per box
☒ -46301 23 G, 6 pcs. per box



Single-Use High Speed Vitreous Cutter
Mach 2 UNO Colorline
double-blade, pneumatic,
with Geuder safety-connector, sterile

☒ -42002 20 G, 1 pc. per box
☒ -46302 23 G, 6 pcs. per box
☒ -42502 25 G, 1 pc. per box
☒ -46502 25 G, 6 pcs. per box



Single-Use High Speed Vitreous Cutter
Mach 4 UNO Colorline
quadruple-blade, pneumatic,
with Geuder safety-connector, sterile

☒ -46304 23 G, 6 pcs. per box
☒ -46504 25 G, 6 pcs. per box




Geuder® Single Use Instruments
Trocar Systems

HATTENBACH - NIKOLIC
Single-Use Trocar System One-Step
UNO Colorline
3 ports for one-step technique,
incl. infusion tubing, sterile

☒ -46311 23 G, 6 pcs. per box
☒ -46511 25 G, 6 pcs. per box








trocar head MVR lancet






	
SZURMAN Single-Use DMEK Cartridge for Descemet Membrane Endothelial Keratoplasty incl. tube connection for loading the cartridge, sterile	
GS-38630	incision 2.8 to 3.0 mm, 14 Gauge/ 2.0 mm
GS-38635	incision 2.4 to 2.75 mm, 16 Gauge/ 1.6 mm
	
SZURMAN DMEK Rapid Transportation Cartridge for transportation of prepared Descemet membrane for use in Descemet Membrane Endothelial Keratoplasty (DMEK) incl. cartridge, 2 plugs, transportation holder for cartridge and connection for loading the cartridge 6 pcs. per box, sterile	
GS-38626	
	
Single-Use Tubing with Luer Connector e.g. for rinsing and staining the Descemet membrane 6 pcs. per box, sterile	
GS-38611	

	
Single-Use Iris Retractor for mechanical dilation of the pupil, flexible 6 iris retractors per set, 5 sets per box, sterile	
GS-34524	
GS-34525	
	
Lacrimal Intubation Tubing silicone tubing 35 cm outside diameter 0.6 mm, inside diameter 0.3 mm 10 pcs. per box, sterile	
GS-34810	




Geuder® Single Use Instruments
Cannulas

	
Single-Use Lacrimal Cannula curved 26 Gauge/ 0.45 x 30 mm 10 pcs. per box, sterile	
GS-34290	
	
ATKINSON Single-Use Retrobulbar Cannula 10 pcs. per box, sterile 25 Gauge, 0.5 x 40 mm 23 Gauge, 0.64 x 38 mm	
GS-34210	
GS-34215	
	
Single-Use Peribulbar Cannula 25 Gauge, 0.5 x 32 mm 10 pcs. per box, sterile	
GS-34220	
	
Single-Use Cystotome Cannula vaulted, 45° 25 Gauge/ 0.5 mm 10 pcs. per box, sterile 12.5 mm long, side cutting 12.5 mm long	
GS-34225	
GS-34230	
	
LEWICKY Single-Use A. C. Maintainer 20 Gauge/ 0.9 mm, tip 5 mm 3 pcs. per box, sterile	
GS-34261	

Geuder® 15 Single-Use Instruments
Cannulas

⌘-34235 ⌘-34240	RVCROFT Single-Use A. C. Cannula angled at 7 mm, 45° 10 pcs. per box, sterile 20 Gauge, 0.9 x 5 + 17 mm 25 Gauge, 0.5 x 25 mm	  
⌘-34245	SAUTTER Single-Use Hydrodissection Cannula vertically flattened 27 Gauge/ 0.4 x 22 mm 10 pcs. per box, sterile	  
⌘-34270	TUEBINGEN MODEL Single-Use Hydrodissection Cannula 30 Gauge/ 0.3 x 19.5 mm 10 pcs. per box, sterile	  
⌘-34271	FRANKFURT MODEL Single-Use Hydrodissection Cannula 27 Gauge/ 0.4 mm x 15 mm 10 pcs. per box, sterile	  
⌘-34500	SAAD Single-Use Unfolding Cannula for unfolding Descemet membrane in DMEK 1 front port, 4 lateral ports 27 Gauge/ 0.4 mm 10 pcs. per box, sterile	  
⌘-34501	Single-Use DALK Cannula for big-bubble technique of Deep Anterior Lamellar Keratoplasty (DALK) 27 Gauge/ 0.4 mm, curved, flattened tip 0.4 mm posterior port 10 pcs. per box, sterile	  
⌘-34265	Single-Use Irrigation Vectis Cannula double port 25 Gauge/ 0.5 mm 10 pcs. per box, sterile	  

Geuder® 15 Single-Use Instruments
Cannulas

⌘-34280	SIMCOE Single-Use I/A Cannula 0.4 mm diameter aspiration port front opening irrigation port 23 Gauge/ 0.64 mm 10 pcs. per box, sterile	  
⌘-80975	Single-Use Injection Cannula for EasyGas® 30 Gauge/ 0.3 x 12 mm 100 pcs. per box, sterile	  
⌘-34400	Single-Use Bimanual I/A Handpiece Set consisting of ■ irrigation handpiece, irrigation ports 2 x 0.6 mm ■ aspiration handpiece, aspiration port 0.4 mm textured tip 21 Gauge/ 0.8 mm 5 pcs. per box, sterile	   irrigation aspiration  
⌘-34289	Single-Use Backflush Handpiece with silicone chamber, Luer-Lock 10 pcs. per box, sterile	 
⌘-34806	Single-Use Silicone Chamber for backflush handpieces ⌘-37000, ⌘-37002, ⌘-31805, 10 pcs. per box, sterile	 

Geuder® 15 Single-Use Instruments

Cannulas

34275	Single-Use Backflush Cannula inside diameter 25 Gauge/ 0.5 mm outside diameter 21 Gauge/ 0.8 mm 35 mm long 10 pcs. per box, sterile			
34291 34296 34299	Single-Use Backflush Cannula blunt tip, 32 mm long 10 pcs. per box, sterile 20 Gauge/ 0.9 mm 23 Gauge/ 0.6 mm 25 Gauge/ 0.5 mm			
34293 34294 34297	Single-Use Backflush Cannula with silicone brush 5 pcs. per box, sterile 20 Gauge/ 0.9 mm 23 Gauge/ 0.6 mm 25 Gauge/ 0.5 mm			
34292 34295 34298	Single-Use Backflush Cannula with silicone tip 5 pcs. per box, sterile 20 Gauge/ 0.9 mm 23 Gauge/ 0.6 mm 25 Gauge/ 0.5 mm			
34285	Single-Use PFCL Cannula for injection of heavy fluids dual bore, coaxial 23 Gauge/ 0.64 mm 10 pcs. per box, sterile			
34493 34494 34495 34496	Single-Use VFI Cannula for Oil I/A for silicone oil injection/aspiration polyimide tip 6 mm, thin-walled 10 pcs. per box, sterile 20 Gauge/ 0.9 mm 23 Gauge/ 0.6 mm 25 Gauge/ 0.5 mm 27 Gauge/ 0.4 mm			

Geuder® 15 Single-Use Instruments

Cannulas

34497 34498	Single-Use Cannula for silicone oil injection 10 pcs. per box, sterile 20 Gauge/ 0.9 mm x 8 mm 23 Gauge/ 0.6 mm x 8 mm			
33468 33469	Single-Use Infusion Cannula with 25 cm silicone tube Luer-Lock plastic adapter 20 Gauge/ 0.9 mm 5 pcs. per box, sterile beveled tip, 4.0 mm beveled tip, 6.0 mm			
33488 33489	Single-Use Injection Cannula for viscous fluids with 1 metal sleeve, Luer-Lock plastic adapter and 25 cm PVC tube 20 Gauge/ 0.9 mm 5 pcs. per box, sterile beveled tip, 4.0 mm beveled tip, 6.0 mm			

Geuder® 15 Single-Use Instruments

Injection Syringes

31891	Single-Use Syringe Luer-Lock, 10 ml, sterile		
38610	Single-Use Syringe BD Emerald™ Luer slip, 5 ml, sterile 100 pcs. per box		

⊠-31520

Single-Use Mini Filter
0.1 µm, sterile

⊗



⊠-23991

Single-Use Three-Way Stop Cock
2x Luer-Lock female,
1x Luer-Lock male plastic,
sterile

⊗




⊠-30146

Single-Use Air Injection Tube
for megaTRON S3/ S4 HPS
incl. filter 0.1 µm and 2-way-stopcock
length 1.9 m
6 pcs. per box, sterile

⊗






⊠-28766

Single-Use Oil Injection System
to inject silicone oil pneumatically,
with protective cover for glass syringe,
pressure tube fits megaTRON® S3/ S4 HPS
and Pentasys® 2, sterile


⊗



⊠-28767

Single-Use Oil Injection System
to inject silicone oil pneumatically,
with protective cover for glass syringe,
pressure tube fits megaTRON® and Accurus®, sterile

⊗



⊠-28768


Single-Use Oil Injection System
to inject silicone oil pneumatically,
with protective cover for glass syringe,
pressure tube fits Millennium®, sterile

⊗

Geuder®

15 Single-Use Instruments







Cautery, Scalpel Handle, Absorbent Swabs, Marking Pen

⌘-19828	Single-Use Low Temperature Battery Cautery fine tip (1.5 mm), tip length 28 mm, max. temp. 600°C 10 pcs. per box, sterile	
⌘-19829	Single-Use High Temperature Battery Cautery fine tip, tip length 28 mm, max. temp. 1,200°C 10 pcs. per box, sterile	
⌘-32700	Scalpel Handle	
		
⌘-18127	Absorbent Swabs 17 mm long box of 20 pouches 10 swabs per pouch sterile	
⌘-33786	Marking Pen gentian violet dual tip, fine and broad to ink corneal instruments 20 pcs. per box, sterile	

Geuder®

15 Single-Use Instruments

Silicone Tires

⌘-34602	276 Style Silicone Tire, Asymmetrical 7.0 mm wide 2.5 mm groove 5 pcs. per box, sterile	 
⌘-34606	280 Style Silicone Tire, Asymmetrical 10.0 mm wide 2.5 mm groove 5 pcs. per box, sterile	 
⌘-34608	286 Style Silicone Tire, Convex 6.0 mm wide 2.5 mm groove 5 pcs. per box, sterile	 
⌘-34610	287 Style Silicone Tire, Convex 7.0 mm wide 2.5 mm groove 5 pcs. per box, sterile	 
⌘-34616	277 Style Silicone Tire, Concave 7.0 mm wide 2.5 mm groove 5 pcs. per box, sterile	 
⌘-34618	279 Style Silicone Tire, Concave 9.0 mm wide 2.5 mm groove 5 pcs. per box, sterile	 

Geuder® 15 Single-Use Instruments


Silicone Circling Bands


34620	40 Style Silicone Circling Band W 2.0 x H 0.75 x L 125 mm 5 pcs. per box, sterile		
	suitable sleeves		
	34642 (round)		
	34645 (oval)		

34622	240 Style Silicone Circling Band W 2.5 x H 0.6 x L 125 mm 5 pcs. per box, sterile		
	suitable sleeves		
	34645 (oval)		
	34646 (round)		

Geuder® 15 Single-Use Instruments

Silicone Strips


34624	41 Style Silicone Strip W 3.5 x H 0.75 x L 125 mm 5 pcs. per box, sterile		
	suitable sleeves		
	34647 (oval)		
	34644 (round)		


34626	42 Style Silicone Strip W 4.0 x H 1.25 x L 125 mm 5 pcs. per box, sterile		
	suitable sleeves		
	34647 (oval)		
	34644 (round)		

34628	4050 Style Silicone Strip W 5.0 x H 1.0 x L 125 mm 5 pcs. per box, sterile		
	suitable sleeves		
	34647 (oval)		
	34644 (round)		

Geuder® 15 Single-Use Instruments


Silicone Strips

34636	219 Style Grooved Silicone Strip W 4.5 x H 1.25 x L 100 mm 2.5 mm groove 5 pcs. per box, sterile		

34638	220 Style Grooved Silicone Strip W 6.0 x H 1.5 x L 100 mm 2.5 mm groove 5 pcs. per box, sterile		

Geuder® 15 Single-Use Instruments

Silicone Sleeves


34646	270 Style Silicone Sleeve diameter 0.76 mm inside diameter 1.65 mm outside 30 mm long 5 pcs. per box, sterile		

34642	70 Style Silicone Sleeve diameter 1.0 mm inside diameter 2.1 mm outside 30 mm long 5 pcs. per box, sterile		

34644	72 Style Silicone Sleeve diameter 1.5 mm inside diameter 2.4 mm outside 30 mm long 5 pcs. per box, sterile		




Geuder® 15 Single-Use Instruments

Silicone Sleeves

34645	3083 Style Oval Silicone Sleeve diameter 1.5 mm inside diameter 2.5 mm outside 5 mm long 5 pcs. per box, sterile	
34647	3084 Style Oval Silicone Sleeve diameter 1.5 mm inside diameter 3.75 mm outside 5 mm long 5 pcs. per box, sterile	






Geuder® 15 Single-Use Instruments





Silicone Sponge Implants

34656	503 Style Round Silicone Sponge Implant diameter 3.0 mm 80 mm long 5 pcs. per box, sterile	
34658	504 Style Round Silicone Sponge Implant diameter 4.0 mm 80 mm long 5 pcs. per box, sterile	
34660	505 Style Round Silicone Sponge Implant diameter 5.0 mm 80 mm long 5 pcs. per box, sterile	

Geuder® 15 Single-Use Instruments

Silicone Sponge, Implants

34664	501 Style Oval Silicone Sponge Implant W 4.0 x H 2.5 x L 80 mm 5 pcs. per box, sterile	
34666	506 Style Oval Silicone Sponge Implant W 5.0 x H 3.0 x L 100 mm 5 pcs. per box, sterile	
34668	507 Style Oval Silicone Sponge Implant W 7.5 x H 5.5 x L 80 mm 5 pcs. per box, sterile	
34672	506G Style Grooved Silicone Sponge Implant 2.5 mm groove W 5.0 x H 3.0 x L 80 mm 5 pcs. per box, sterile	
34678	509G Style Grooved Silicone Sponge Implant 2.5 mm groove W 7.5 x H 3.5 x L 80 mm 5 pcs. per box, sterile	

34680	516G Style Grooved Silicone Sponge Implant 2.5 mm groove W 6.4 x H 2.3 x L 80 mm 5 pcs. per box, sterile	
34684	510 Style Semi-Circular Silicone Sponge Implant W 5.0 x H 2.5 x L 80 mm 5 pcs. per box, sterile	
34686	511 Style Semi-Circular Silicone Sponge Implant W 7.5 x H 2.75 x L 80 mm 5 pcs. per box, sterile	
34688	515 Style Semi-Circular Silicone Sponge Implant W 5.0 x H 3.33 x L 80 mm 5 pcs. per box, sterile	

34711	Single-Use Set for Intravitreal Injection (IVOM) consisting of - measuring and fixation instrument, 3.5/ 4.0 mm - eye speculum - irrigation cannula (0.8 mm x 60 mm) - cannula 0.9 x 70 mm - syringe 1 ml - syringe 10 ml - eye drape with pre-cut adhesive incise-film 50 x 50 cm opening 7 x 10 cm - absorbent swab (3x) - gauze compress 10 x 10 cm (2x) - protective cover 60 x 60 cm 30 pcs. per box, sterile	
-------	--	---