

Geuder 16
Miscellaneous

Sonstiges



- Capsular Tension Ring Injectors 199
- Lacrimal, DCR Instruments 199
- Silicone Blocks 202
- Vitrectomy Lenses 203
- Eye Shield, Magnet 204
- Diathermy Instruments 205

Geuder 16 Miscellaneous
Capsular Tension Ring Injectors

G-32955

Injector
for capsular tension ring
shaft diameter 1.2 mm
for 1.9 mm incisions
shaft length 30 mm
overall length 156 mm



G-32960

Injector
for MORCHER capsular tension ring
shaft diameter 1.2 mm
for 1.9 mm incisions
shaft length 30 mm
overall length 156 mm



Teflon coated for
smooth gliding



Geuder 16 Miscellaneous
Lacrimal, DCR Instruments

G-31810

FRIEDMANN
Rongeur
overall length 140 mm



Scale 2:3

G-32355

Nasal Speculum
overall length 150 mm



Scale 2:3

G-33125

MEAD
Hammer
300 g, lead filled
head diameter 27 mm
overall length 190 mm



Scale 2:3

G-18630

HEIDELBERG MODEL
Elevatorium
overall length 143 mm



G-18620

HEIDELBERG MODEL
Raspatory
overall length 155 mm



G-18240

WEST
Bone Chisel
for lacrimal sac
flat, 5.0 mm wide
overall length 160 mm



G-18245

WEST
Bone Gouge
for lacrimal sac
concave, 5.0 mm wide
overall length 162 mm



G-33652

LASEK Trephine
with crosshairs
9.0 mm
overall length 100 mm




G-33655
G-33657

LOHMANN
LASEK Funnel
with crosshairs
overall length 100 mm
8.5 mm
9.5 mm



G-32313

WRIGHT
Fascia Needle
overall length 115 mm



Geuder 16 Miscellaneous
Silicone Blocks

GS-32671	ÖRK LINZ MODEL Silicone Block for DMEK for donor cornea white-colored outside diameter 30 mm height 19 mm only for single use, not sterile	
GS-32672	Silicone Block for donor cornea blue-colored outside diameter 30 mm height 19 mm only for single use, not sterile	
GS-32674	Silicone Block with center point for donor cornea blue-colored outside diameter 20 mm height 10 mm only for single use, not sterile	
GS-38620	SZURMAN DMEK Graft Retainer for preparation of DMEK graft white-colored outside diameter 35 mm height 10 mm Teflon, for multiple use	

Geuder 16 Miscellaneous
Vitrectomy Lenses

GS-36000	KILP Lens 60 D flat, field of vision 20°, silicone ring central posterior pole and central vitreous body	
GS-36010	LANDERS Lens 93 D, biconcave with silicone ring clear view of fundus in air-filled phakic eye	
GS-36015	MACHEMER Lens with silicone ring to magnify structures in vitreous cavity or on retinal surface	
GS-36020	PEYMAN Lens wide field of vision 30°, with silicone ring for central fundus and central vitreous body	
GS-36025	TOLENTINO Lens prism 20° with silicone ring posterior peripheral fundus and posterior peripheral vitreous body	
GS-36030	TOLENTINO Lens prism 30° with silicone ring peripheral viewing beyond the equator	
GS-36045	Silicone Replacement Ring for KILP Lens GS-36000	
GS-36050	Silicone Replacement Ring for lenses from GS-36010 to GS-36035	
GS-36065	Lens Handle Ring for lenses without silicone ring to suture on to sclera for GS-36000	

GS-31105	Eye Shield plastic, transparent for use on both eyes 100 pcs. per bag	
GS-31900	KUHN Intraocular Magnet 30 mm shaft, 1.9 mm outside diameter 4 mm magnetic tip overall length 90 mm	

GS-31215	Bipolar Forceps straight, 100 mm (sterilizable in autoclave only)	
GS-31216	Bipolar Forceps angled, 100 mm (sterilizable in autoclave only)	
GS-31218 GS-31219	Bipolar Handpiece diameter 1.1 mm diameter 0.9 mm	
GS-25030 GS-25032 GS-25034	Intraocular Diathermy Probe 20 Gauge/0.9 mm 23 Gauge/0.64 mm 25 Gauge/0.5 mm	
GS-28220	Diathermy Connection Cable	