

17 Sterilization Containers and Adapters

Sterilisationscontainer und Reinigungsadapter

- microStop® Containers and Accessories 207
- Small Sterilization Cases 209
- Sterisafe® Containers and Accessories 209
- Silicone Strips, Silicone Mat, Seals 210
- Small Cases 211
- Rinsing and Cleaning Adapters 213



Geuder 17 Sterilization Containers and Adapters microStop® Containers and Accessories

- 39783
- 39786

microStop® Sterilization Tray Mini
29 mm, low
L 277 x W 172 x H 29 mm
without lid
with lid




- 39785
- 39784

microStop® Sterilization Tray Mini
38 mm, high
L 277 x W 172 x H 38 mm
without lid
with lid



- 39787

microStop® Sterilization Tray Lid Mini
L 277 x W 172 mm
for 39783 and 39785



- 39782
- 39789

microStop® Retainer Element
for microStop® tray mini
nap style, 120 mm
separable into 6 single elements
29 mm, low
38 mm, high



- 39788

microStop® Clip for Retainer Elements
for 39782 and 39789
10 pcs. per bag



☒-39741 microStop® Label
for lid, customized labeling



☒-39746 microStop® ID Label
for front, customized labeling




☒-39753 color red
☒-39754 color blue
☒-39755 color green

microStop® Logistic Frame
for lid
to use with labels




☒-39762 microStop® Seal
with indicator
1000 pcs. per box



☒-39764 microStop® Tag
with indicator
1000 pcs. per box




☒-31649 Indicator
for steam sterilization, slim
1000 pcs. per box



☒-31657 Safety Seal
for steam sterilization cases
100 pcs. per box




☒-33100 Sterilization Case
for 1 vitreoretinal instrument or diamond knife
plastic, with 2 bridges (steam sterilization)
outside dimensions 195 x 60 x 25 mm



☒-33101 Sterilization Case
for 1 vitreoretinal instrument
plastic, with burl mat (steam sterilization)
outside dimensions 195 x 60 x 25 mm



☒-33102 Sterilization Case
for 4 – 5 instruments
plastic, with burl mat (steam sterilization)
outside dimensions 200 x 110 x 28 mm



☒-39527 Sterisafe® Duro A4 Plus
container for steam, gas (EO, F0) and
plasma sterilization, without trays and strips
outside dimensions 290 x 190 x 70 mm
inside dimensions: 265 x 165 x 60 mm



☒-39525

Sterisafe® Duro A4

container for steam, gas (EO, FO)
and plasma sterilization
outside dimensions 290 x 190 x 60 mm
inside dimensions: 265 x 165 x 49 mm



☒-39528

Sterisafe® Duro A6

container for steam, gas (EO, FO)
and plasma sterilization
outside dimensions 190 x 110 x 60 mm
inside dimensions: 165 x 92 x 49 mm



☒-39543

Sterisafe® Duro Wire Basket A4-30

outside dimensions 260 x 163 x 30 mm
inside dimensions: 249 x 157 x 23 mm



☒-39541

Sterisafe® Duro Metal Tray A4-15

outside dimensions 260 x 163 x 15 mm
inside dimensions: 249 x 157 x 11 mm



☒-39542

Sterisafe® Duro Metal Tray A4-30

outside dimensions 260 x 163 x 30 mm
inside dimensions: 249 x 157 x 26 mm



☒-39549

Sterisafe® Duro Wire Basket A6-30


outside dimensions 162 x 90 x 30 mm
inside dimensions: 156 x 80 x 23 mm



☒-39548

Sterisafe® Duro Metal Tray A4-15

outside dimensions 162 x 90 x 15 mm
inside dimensions: 156 x 80 x 11 mm



☒-39561

Sterisafe® Silicone Fixation Strip 75

with spring action
for Sterisafe® Duro A4 and A6
length: 75 mm
5 pcs. per box



☒-39562

Sterisafe® Silicone Fixation Strip 145

with spring action
for Sterisafe® Duro A4 and A6
length: 145 mm
5 pcs. per box



☒-39563

Sterisafe® Silicone Fixation Strip 230

with spring action
for Sterisafe® Duro A4
length: 230 mm
5 pcs. per box



☒-39501
☒-39511

Sterisafe® Filter

A4 10 pcs. per box
A6 10 pcs. per box




☒-39530

Sterisafe® Filter Tag Expiration Date

10 pcs. per box




☒ -39535	Seal/ Label seal with color indicator and blank label 2000 pcs. per box	<table><tr><td>Steri. - Datum/Date</td><td></td></tr><tr><td>Name</td><td></td></tr><tr><td></td><td>Dampf (braun)  Steam (brown)</td></tr></table>	Steri. - Datum/Date		Name			Dampf (braun)  Steam (brown)
Steri. - Datum/Date								
Name								
	Dampf (braun)  Steam (brown)							
☒ -39564	Universal Label Holder for Sterisafe® Duro 2 pcs. per box							
☒ -39540	Security Seal 100 pcs. per box							
☒ -39502 ☒ -39512	Sterisafe® Silicone Burl Mat A4 A6							

⌘-39544	Rinsing Adapter for Wire Baskets including 4 silicone tubings 2x male Luer adapters 2x female Luer-Lock adapters and 2 screws	
⌘-39557	Silicone Tubings for ⌘-39544 Set of 4 pcs. 2x male Luer adapters 2x female Luer-Lock adapters 10 sets per box	
⌘-31028	Cleaning Adapter for ultrasonic handpieces (irrigation connector) Luer female/ male	
⌘-31029	Cleaning Adapter for ultrasonic handpieces (aspiration connector) Luer female/ male	
⌘-31021	Cleaning Adapter for ultrasonic tips	



┐-31013 **Cleaning Adapter**
for HATTENBACH 25 Gauge Hybrid Instruments
(┐-32944, ┐-32947 and ┐-32948)



┐-22636 **Adapter**
for PAPST I/A handpiece
light model



┐-31025 **Cleaning Adapter**
metal
for cleaning I/A handpieces
Luer-Lock female/ female



┐-31026 **Cleaning Adapter**
plastic
for cleaning I/A handpieces
Luer-Lock female/ female



┐-31022 **Cleaning Adapter**
for ┐-37010 to ┐-37080



┐-31031 **Cleaning Adapter**
for Aspiration Handpiece ┐-37100
Luer-Lock female/ silicone tube
2.1 x 1.0 mm



┐-31018 **Cleaning Adapter**
for PFCL Cannula ┐-33055
Luer-Lock female/ silicone tube



┐-31024 **Cleaning Adapter**
for eye specula with aspiration



┐-31023 **Cleaning Adapter**
for lacrimal cannulas with handle
and lateral aspiration port 0.40 mm



┐-31027 **Cleaning Adapter**
for JUENEMANN probes



GS-31034 **Cleaning Adapter**
for foreign body instruments GS-14540, GS-14550
Luer-Lock female

GS-31014 **Cleaning Adapter**
for Cannulas GS-33057
Luer-Lock female



GS-36604 **Cleaning Adapter**
for rotatable and interchangeable Ergoline tips
GS-36902, GS-36915, GS-36925, GS-36930, GS-36955



Please read your instrument's cleaning instructions
or ask your local sales representative
for the right cleaning adapter.