

SINGLE-USE PHACO TIPS

**FOR COAXIAL
PHACOEMULSIFICATION**




- High followability and holdability
- Short ultrasound times with minimal stress due to optimized inner diameter
- Maximized active surface for ultrasound emulsification
- Improved quality of tip surface and shape
- Gentle insertion due to smooth transition to sleeve
- Enhanced sleeve irrigation

SINGLE-USE PHACO KITS

including:

- ultrasonic tip, titanium, 30°
- infusion sleeve
- test anterior chamber
- wrench
- 10 pcs per box, sterile



STEPPED (FLARED) TIPS

G-24299
for 2.2 - 2.35 mm
incision



G-24303
for 2.4 - 2.55 mm
incision



G-24306
for 2.6 - 2.75 mm
incision



G-24307
for 2.8 - 2.95 mm
incision



STRAIGHT TIPS

G-24351
for 2.2 - 2.35 mm
incision



G-24352
for 2.4 - 2.55 mm
incision



G-24353
for 2.6 - 2.75 mm
incision



G-24354
for 2.8 - 2.95 mm
incision



GEUDER AG reserves the right to make changes to technical details in response to recent developments. GEUDER does not assume liability for the accuracy of each individual statement.

Illustrations not drawn to scale.